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Part 2. Usage Instruction

2.1. General provisions.....	43
2.2. Conditions for treatment.....	45
2.3. Course of treatment.....	45

<i>Appendix. Recommended method of use of External Paraorbital Electrode in Ophthalmology.....</i>	<i>47</i>
--	-----------

Warranty Certificate.....	81
---------------------------	----

Acceptance Certificate.....	83
-----------------------------	----

PART 1

Technical Passport

1.1. GENERAL DESCRIPTION

External paraorbital electrode (EPE) is intended for preventive and therapeutic noninvasive impact on biologically active zones around eyes by the method of dynamic electrostimulation.

Dynamic electric neurostimulation facilitates recovery from visual impairments caused by some eye diseases. Additionally, electric neurostimulation has a general regulating effect on the physiological systems of human organism.

The unit is intended for household and clinical use.

The External Paraorbital Electrode should only be used in combination with the DENAS and DiaDENS (DiaDENS-T, -PC, -PCM)* electric stimulation devices.

** Connection of the electrode to some device models with old-style sockets should be made via an adapter (which is not included in this complication). Models of DENAS devices without sockets for external electrode connection require a special support (not included in this complication).*

Compliance with standards: *This medical device is CE marked according to the Directive 93/42 /EEC on medical equipment.*

Attention! *If you connect the electrode to the models of devices produced before 2009, for proper work it is recommended to disable the sleep mode.*

1.2. TECHNICAL SPECIFICATIONS

1.2.1. Weight of the device is no more than 300 g;

1.2.2. External dimensions are no more than 181x76x57 mm;

1.2.3. Length of the electrode cable is not less than 600 mm.

1.3. PACKAGE CONTENTS

EN

The package includes:

Name	Quantity
External Paraorbital Electrode DENAS-Glasses	1
Connecting cable	1
Operation manual	1
Package	1

1.4. SAFETY RULES

Please read the following safety instructions! They are designed to prevent user and property of damage. The manual contains important recommendations on proper use and care of the unit.

Patients with implanted electronic devices (such as heart pacemaker) and patients with idiosyncrasy to electric current are not allowed to use this unit.

Warning: Do not operate the unit while using another electrical medical device. This can cause burns and irreparable damage to the device. Before you begin stimulation, make sure the patient is not connected to any high-frequency electronic device!

The product contains fragile elements. Prevent it from shocks.

The unit is not waterproof. Keep away from liquids and moisture.

All work for repair of the product must be performed by qualified specialists at the company-manufacturer.

Transportation conditions: temperature from -50 to +50 °C, relative humidity from 30 to 93 %, atmospheric pressure from 70 to 106 kPa (from 525 to 795 mmHg).

Storage conditions: temperature from -50 to +40 °C, relative humidity from 30 to 93 %, atmospheric pressure from 70 to 106 kPa (from 525 to 795 mmHg).

Operating conditions temperature from +10 to +35 °C, relative humidity from 30 to 93 %, atmospheric pressure from 70 to 106 kPa (from 525 to 795 mmHg).

Attention! *If the unit was stored at the ambient temperature below 10 °C keep it under normal climatic conditions for at least two hours before use.*

Recycling: the unit packaging materials are not hazardous for the environment; they can be recycled.

1.5. TECHNICAL MAINTENANCE

Maintenance of the unit should include the following:

- visual check of the device;
- cleaning of electrodes (use standard disinfection detergents and a soft cloth to clean electrodes);
- check of the unit's readiness for operation when connected to the DENAS or DiaDENS devices.

1.6. WARRANTY

1.6.1. The manufacturer guarantees the compliance of the device with Technical Conditions, provided the terms of operation, transportation and storage are observed.

1.6.2. The operation lifetime of the unit is 5 years.

Observation of regulations governing the device application can considerably increase the lifetime set by the manufacturer.

1.6.3. This product is warranted for the period of 24 months from the original date of purchase.

1.6.4. The consumer is to check the package contents and perform the visual control of the device in the presence of the retailer. Complaints about the incomplete set and appearance of the device are not accepted after sale.

1.6.5. In case of identification of defects during warranty period, the seller, manufacturer or manufacturer's representative specified by consumer rights legislation. The seller, manufacturer or manufacturer's representative is not responsible for the defaults should they occur after sale of the device in the result of:

- 1) failure on the part of the consumer to comply with the rules of transportation, storage, care and operation, provided by the present manual;
- 2) actions of third party;

3) force-majeure circumstances.

EN

1.6.6. If warranty servicing is required (i.e., in case of some parts missing or malfunction of the device), an application for repair or replacement should be fully and properly filled in by the customer and presented to the manufacturer. The application must include the following information: customer's name, address, telephone number, brief description of the problem, time and conditions of its occurrence.

Warranty repairs must be carried out by the Manufacturer or Manufacturer's Service Centers.

1.7. MANUFACTURER'S ADDRESS

TRONITEK LLC

620146, Russia, Ekaterinburg, Akademika Postovskogo Str., 15.

Telephone: +7 (343) 267-23-30

<http://www.denascorp.ru>, e-mail: corp@denascorp.ru

Official representative in EU countries:

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Usage instruction

2.1. GENERAL PROVISIONS

The majority of eye diseases are evidence of adaptive dysfunction. The course of an eye disease depends on the patient's general well-being and presence of concomitant conditions.

Dynamic electric neurostimulation (DENS) delivered by the External Paraorbital Electrode is regarded as a constituent part of the prophylactic and therapeutic program in ophthalmopathy treatment.

The stimulation of the paraorbital area has a healthy influence on visual functions as well as optimizes general health.

All patients should be supervised by an ophthalmologist for making a diagnosis, adjusting multimodality treatment and monitoring and evaluating the results.

Contraindications to treatment with the EPE are described in Operation Manuals and User's Instructions to the DENAS devices.

2.2. CONDITIONS FOR TREATMENT


The application of the External Paraorbital Electrode doesn't require special conditions. During the session the patient should sit or lie in a comfortable position.

2.3. COURSE OF TREATMENT

2.3.1. Prior to the start of treatment, wipe the surface of the electrode with a soft cloth moistened with a disinfecting liquid, (e.g. 3 % solution of hydrodgen peroxide).

2.3.2. Put on and anchor the unit.

2.3.3. Connect EPE to the DENAS device. If the socket of the device does not match the plug of the electrode, use the adapter.

2.3.4. Turn the device on.

2.3.5. Choose the certain mode and set the individual power (intensity of impact) (see appendix).

2.3.6. Conduct the session. The duration of the procedure with the use of EPE is defined by the patient's age, disease and the intensity of symptoms (see appendix).

2.3.7. Upon termination of the session turn the device off, remove and disconnect EPE from the device.

Attention! *It is recommended to switch the workout of segmental and reflexogenic zones in the schemes which are used to treat chronic eye diseases. Choice and treatment of extra zones by inbuilt electrodes apparatuses in accordance general rules which are described in the Operation Manual of the appropriate unit*.*

*Guidance manual on dynamic electrostimulation. Ekaterinburg, 2011.

Attention! Turn the EPE off during applying this device to another zones.

APPENDIX

Recommended method of use of External Paraorbital Electrode in Ophthalmology

Diagnosis	Device	Modes	Session				Note
			Duration	Stimulation frequency (Hz)	Power	Quantity	
Asthenopia (visual fatigue against the continuous and intensive stress at reading, writing, computer work, etc.)	DENAS	'Therapy'	3–5 minutes	77 Hz	minimal	1 session per day (in the evening), the course of 10–14 days	Repeated courses are conducted against the intensive visual stresses when necessary or 3–4 courses per year
	DiaDENS	'Therapy'	3–5 minutes	10, 20 Hz or 77 Hz			

Diagnosis	Device	Modes	Session			Note	
			Duration	Stimulation frequency (Hz)	Power		Quantity
Presbyopia (age-related far-sightedness)	DENAS	'Therapy'	3–5 minutes	77 Hz (Th 1)	mini-mal	1 session per day, the course of 10–14 days	Repeated course is conducted in 20–30 days. Take a pause of 2–3 months and repeat the courses cycle, up to 3–4 courses per year
	DiaDENS	'Therapy'	3–5 minutes	10, 20 Hz or 77 Hz			

Diagnosis	Device	Modes	Session				Note
			Duration	Stimulation frequency (Hz)	Power	Quantity	
Spasm of accommodation Myopia (short-sightedness) Hypermetropia Astigmatism Amblyopia (weak eyesight)	DENAS	'Therapy'	children up to 7 years – 3 min. 8–12 years – 5 min., over 12 years 5–7 min.	77 Hz (Th 1)	minimal	1 session per day, the course of 10–14 days	Repeated courses 3–6 times per year
	DiaDENS	'Therapy'		20 Hz or 77 Hz			

Diagnosis	Device	Modes	Session				Note
			Duration	Stimulation frequency (Hz)	Power	Quantity	
<i>Cataract</i>	DENAS	'Therapy'	3–5 min.	77 Hz (Th 1)	minimal	1 session per day, the course of 7–10 days	Conduct 4–6 courses with an interval of 10–20 days
	DiaDENS	'Therapy'	3–5 min.	20 Hz or 77 Hz			
<i>Glaucoma</i>	DENAS	'Therapy'	3–5 min.	77 Hz (Th 1)	minimal	1 session per day, the course of 10 days	Repeated courses when necessary, depending on the dynamics of internal eye pressure
	DiaDENS	'Therapy'	3–5 min.	77 Hz			

Diagnosis	Device	Modes	Session				Note
			Duration	Stimulation frequency (Hz)	Power	Quantity	
<i>Inflammatory eye diseases</i> (at acute disease and exacerbation of chronic process)	DENAS	'Therapy'	Sessions of 5–10 min. with the pause of 1–1.5 hours till the symptoms of inflammation decrease	77 Hz (Th 1)	mini-mal	After the improvement conduct 1 session per day until the recovery	Before and after each session it is necessary to treat EPE with a disinfectant solution
	DiaDENS	'Therapy'		77 Hz			

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Signature of the official of the company responsible for the acceptance after repairs

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I have checked the unit and have no complaints about its appearance and package set

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External Paraorbital Electrode «DENAS-Glasses» is found to be fit for service.

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I hereby confirm that I agree with and accept warranty terms and conditions.

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Carefully examine the electrode before the purchase! Defects of the body (scratches, cracks, chips) are not guarantee events. EPE with such defects is not subject to replacement, repairs or return by warranty.

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